

## Anti-Recombinant-alpha Tubulin antibody [SM5069] (STJA0005069)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/IF/ICC      |
| Host / Source | Mouse              |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                          |
|---------------------|--------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                               |
| Conjugation         | Unconjugated                                                             |
| Purification        | Affinity purification                                                    |
| Dilution Range      | WB (H, M, R) 1:3000-1:6000<br>IHC/IF (H, M, R) 1:1000-1:2000             |
| Formulation         | PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.           |
| Isotype             | IgG1                                                                     |
| Storage Instruction | Store at -20C for up to one year, and avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION


|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">7277</a>        |
| Gene Symbol | <a href="#">TUBA4A</a>      |
| UniProt ID  | <a href="#">TBA4A_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraformaldehyde-fixed HeLa using alpha Tubulin at dilution of 1:1000



Immunofluorescent analysis of paraformaldehyde-fixed rat colon using alpha Tubulin at dilution of 1:2000