

Anti-Lamin A + Lamin C antibody [SM5079] (STJA0005079)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:1000-1:2000 IHC/IF (H, R) 1:1500-1:3000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG2b
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION





Gene ID	16905
Gene Symbol	Lmna
UniProt ID	LMNA_MOUSE

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of Lamin A/C at dilution of 1:2000



Immunohistochemistry analysis of paraffin-embedded mouse liver using Lamin A/C at dilution of 1:3000



Immunohistochemistry analysis of paraffin-embedded rat liver using Lamin A/C at dilution of 1:3000