

Anti-Calbindin antibody [SM5086] (STJA0005086)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Mouse/Rat

PRODUCT PROPERTIES



Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (M, R) 1:500-1:1000 IHC/IF (H, M, R) 1:500-1:1000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	12307
Gene Symbol	Calb1
UniProt ID	CALB1_MOUSE

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Immunohistochemistry of paraffin embedded human kidney using calb1 at dilution of 1:500




Immunohistochemistry of paraffin embedded human renal cancer using calb1 at dilution of 1:500




Immunohistochemistry of paraffin embedded mouse brain using calb1 at dilution of 1:500



Immunohistochemistry of paraffin embedded mouse kidney using calb1 at dilution of 1:500



Immunohistochemistry of paraffin embedded rat brain using calb1 at dilution of 1:500



Immunohistochemistry of paraffin embedded rat kidney using calb1 at dilution of 1:500