

Anti-FOXO1 antibody [SM5088] (STJA0005088)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H, M, R) 1:500-1:2000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG2b
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	56458
Gene Symbol	Foxo1
UniProt ID	FOXO1_MOUSE

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




Immunohistochemistry of paraffin embedded human colon using FOXO1 at dilution of 1:1000




Immunohistochemistry of paraffin embedded human esophagus cancer using FOXO1 at dilution of 1:1000




Immunohistochemistry of paraffin embedded mouse heart using FOXO1 at dilution of 1:1000




Immunohistochemistry of paraffin embedded mouse lung using FOXO1 at dilution of 1:1000




Immunohistochemistry of paraffin embedded mouse testis using FOXO1 at dilution of 1:1000




Immunohistochemistry of paraffin embedded rat liver using FOXO1 at dilution of 1:1000



Immunohistochemistry of paraffin embedded rat pancreas using FOXO1 at dilution of 1:1000



Immunohistochemistry of paraffin embedded rat spleen using FOXO1 at dilution of 1:1000



Immunohistochemistry of paraffin embedded rat testis using FOXO1 at dilution of 1:1000