

Anti-HDAC1 antibody [SM5093] (STJA0005093)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:500-1:1000 IHC/IF (H, M, R) 1:1000-1:2000/1:1000-1:2000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	433759
Gene Symbol	Hdac1
UniProt ID	HDAC1_MOUSE

ADDITIONAL INFORMATION


Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---




Western blot analysis of HDAC1 at dilution of 1:1000




Immunohistochemistry analysis of paraffin embedded human lung cancer using HDAC1 at dilution of 1:1000




Immunohistochemistry analysis of paraffin embedded human liver cancer using HDAC1 at dilution of 1:1000




Immunohistochemistry analysis of paraffin embedded human colon cancer using HDAC1 at dilution of 1:1000




Immunohistochemistry analysis of paraffin embedded human ovary cancer using HDAC1 at dilution of 1:1000




Immunohistochemistry analysis of paraffin embedded mouse heart using HDAC1 at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded mouse testis using HDAC1 at dilution of 1:1000



Immunohistochemistry analysis of paraffin embedded rat testis using HDAC1 at dilution of 1:1000




Immunohistochemistry analysis of paraffin embedded mouse brain using HDAC1 at dilution of 1:2000




Immunohistochemistry analysis of paraffin embedded rat brain using HDAC1 at dilution of 1:1000




Immunohistochemistry analysis of paraffin embedded mouse kidney using HDAC1 at dilution of 1:1000




Immunofluorescent analysis of paraformaldehyde-fixed human small intestine cancer using HDAC1 at dilution of 1:1000




Immunofluorescent analysis of paraformaldehyde-fixed rat brain using HDAC1 at dilution of 1:1000




Immunofluorescent analysis of paraformaldehyde-fixed rat colon using HDAC1 at dilution of 1:1000




Immunofluorescent analysis of paraformaldehyde-fixed rat brain using HDAC1 at dilution of 1:1000



Immunofluorescent analysis of paraformaldehyde-fixed rat kidney using HDAC1 at dilution of 1:1000



Immunofluorescent analysis of paraformaldehyde-fixed rat lung using HDAC1 at dilution of 1:1000



Immunofluorescent analysis of paraformaldehyde-fixed rat testis using HDAC1 at dilution of 1:1000