

Anti-Chk2 antibody [SM5109] (STJA0005109)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H) 1:300-1:1000/1:100-1:300
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	50883
Gene Symbol	Chk2
UniProt ID	CHK2_MOUSE

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry of paraffin embedded human colon cancer using CHK2 at dilution of 1:500



Immunohistochemistry of paraffin embedded human colon using CHK2 at dilution of 1:500



Immunohistochemistry of paraffin embedded human stomach cancer using CHK2 at dilution of 1:500