

Anti-eGFP antibody [SM5111] (STJA0005111)

GENERAL INFORMATION



Product Type	Primary antibodies
Applications	WB/ICC/IF
Host / Source	Mouse
Reactivity	eGFP

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:500-1:1000 IHC/IF (Species independent) 1:2000-1:4000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunocytochemistry analysis of paraffin embedded CHO-S using eGFP at dilution of 1:1500 Blue: DAPI Red: Anti-eGFP