

## Anti-Desmin antibody [SM5174] (STJA0005174)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES





Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (M, R) 1:1000-1:2000 IHC/IF (M, R) 1:200-1:1000/1:200-1:1000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">13346</a>
Gene Symbol	<a href="#">Des</a>
UniProt ID	<a href="#">DESM_MOUSE</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry of paraffin embedded rat heart+LPS using Desmin at dilution of 1:500



Immunohistochemistry of paraffin embedded rat skeletal muscle using Desmin at dilution of 1:500