

Anti-FITC antibody [SM5205] (STJA0005205)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	FITC

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	None
Dilution Range	Elisa (Non) 1:5000-1:10000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG1k
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

