




Anti-smooth muscle Myosin heavy chain 11 antibody [SM5215] (STJA0005215)


GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC/IF
	Host / Source	Mouse
	Reactivity	Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB (M, R) 1:500-1:1000 IHC/IF (H, M, R) 1:1000-1:4000/1:500-1:1000
	Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
	Isotype	IgG
	Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.
TARGET INFORMATION		
	Gene Symbol	Myh11 Uniprot ID= MYH11_MOUSE "
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Immunohistochemistry analysis of paraffin embedded human cancer of biliary duct using MYH11 at dilution of 1:1000




Immunohistochemistry analysis of paraffin embedded human esophagus cancer using MYH11 at dilution of 1:1000




Immunohistochemistry analysis of paraffin embedded mouse esophagus using MYH11 at dilution of 1:1000




Immunofluorescent analysis of paraffin embedded human ovarian cancer using MYH11 at dilution of 1:500




Immunofluorescent analysis of paraffin embedded mouse uterus using MYH11 at dilution of 1:500



Immunofluorescent analysis of paraffin embedded human carcinoma of urinary bladder using MYH11 at dilution of 1:500



Immunofluorescent analysis of paraffin embedded mouse breast using MYH11 at dilution of 1:500



Immunofluorescent analysis of paraffin embedded rat bladder using MYH11 at dilution of 1:500