

Anti-Phospho-ABL1-Tyr-412 antibody (STJA0005240)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 IHC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25
Gene Symbol	ABL1
UniProt ID	ABL1_HUMAN
Specificity	This antibody was cross-adsorbed to phospho-tyrosine coupled to agarose and to phospho-c-Abl (Tyr-245) peptide before affinity purification using phospho-c-Abl (Tyr-412) peptide. On SDS-PAGE immunoblots of K-562 treated with pervanadate the antibody

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---