

Anti-ACTN4 antibody (N-Term) (STJA0005245)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81
Gene Symbol	ACTN4
UniProt ID	ACTN4_HUMAN
Immunogen Region	N-Term
Specificity	The antibody detects a 100 kDa* protein corresponding to the molecular mass of Alpha-actinin 4 on SDS-PAGE immunoblots of human A431 cells and rabbit spleen fibroblasts.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---