

Anti-ANXA5 antibody [M017] (STJA0005259)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:1000 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	308
Gene Symbol	ANXA5
UniProt ID	ANXA5_HUMAN
Specificity	Clone M017 detects a 35 kDa* protein corresponding to the molecular mass of Annexin A5 on SDS-PAGE immunoblots of human A431, MDA-MB-231, MeWo, and PC3 cells, as well as human skin, breast, lung, and brain. The antibody works for ELISA, western blot,

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---