

Anti-APC antibody (central Region) [KT44] (STJA0005261)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rat
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:250 ICC 1:50
Formulation	PBS + 0.05% NaN ₃
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene Symbol	Apc
UniProt ID	APC_MOUSE
Immunogen Region	central Region
Specificity	This antibody detects full length APC at 300 kDa* and proteolytic fragments of APC at lower molecular weights. The antibody detects APC or its fragments in SW480, CaCO ₂ , DLD-1, and Hct116 cells. In immunocytochemistry, the antibody detects in APC in

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------