

Anti-ACTR3 antibody (C-Term) (STJA0005268)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IHC
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat/Xenopus

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 IHC 1:100 ICC 1:250
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10096
Gene Symbol	ACTR3
UniProt ID	ARP3_HUMAN
Immunogen Region	C-Term
Specificity	The antibody detects a 47 kDa* doublet corresponding to the molecular mass of Arp3 on SDS-PAGE immunoblots of human A431, Jurkat, and HeLa cells, and in rat PC12 cells.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---