

Anti-Fbxo32 antibody (N-Term) (STJA0005276)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 IHC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	67731
Gene Symbol	Fbxo32
UniProt ID	FBX32_MOUSE
Immunogen Region	N-Term
Specificity	This antibody detects a 41 kDa* protein corresponding to the apparent molecular mass of Atrogin on SDS-PAGE immunoblots of mouse muscle tissue lysates. Similar results were seen in C2C12 cell lysate and the band observed in western blot is specific

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---