

Anti-Phospho-CDH2-Tyr-820 antibody (STJA0005287)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1000
Gene Symbol	CDH2
UniProt ID	CADH2_HUMAN
Specificity	This antibody was cross-adsorbed to phospho-N-Cadherin (Tyr-860) and unphosphorylated N-cadherin (Tyr-820) peptides before affinity purification using phospho-N-cadherin (Tyr-820) peptide. The purified antibody detects a 130 kDa* band corresponding t

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---