

## Anti-CDH3 antibody (N-Term) [M596] (STJA0005291)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000 ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1001</a>
Gene Symbol	<a href="#">CDH3</a>
UniProt ID	<a href="#">CADH3_HUMAN</a>
Immunogen Region	N-Term
Specificity	The P-cadherin antibody detects a 120 kDa* protein corresponding to endogenous P-Cadherin expressed in human A431 cells.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---