

## Anti-CDH5 antibody (C-Term) (STJA0005292)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1003</a>
Gene Symbol	<a href="#">CDH5</a>
UniProt ID	<a href="#">CADH5_HUMAN</a>
Immunogen Region	C-Term
Specificity	The antibody detects a 140 kDa* band corresponding to VE-cadherin in western blots of human endothelial cells, and this reactivity can be specifically blocked using VE-cadherin peptide (CX2235).

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---