

Anti-Phospho-CDH5-Tyr-685 antibody (STJA0005293)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1003
Gene Symbol	CDH5
UniProt ID	CADH5_HUMAN
Specificity	This antibody was cross-adsorbed to an unrelated phospho-tyrosine peptide and unphosphorylated VE-cadherin (Tyr-685) peptide before affinity purification using phospho-VE-cadherin (Tyr-685) peptide. The purified antibody detects a 140 kDa* band corre

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---