

Anti-CANX antibody (N-Term) [M437] (STJA0005294)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000 ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	821
Gene Symbol	CANX
UniProt ID	CALX_HUMAN
Immunogen Region	N-Term
Specificity	This antibody detects a 90 kDa* protein corresponding to the apparent molecular mass of Calnexin on SDS-PAGE immunoblots of human HeLa and rat A7r5 cells. In immunocytochemistry, anti-Calnexin specifically stains endoplasmic reticulum and related ves

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---