

Anti-Phospho-JUP-Tyr-550 antibody (STJA0005309)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/ICC
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Conjugation	Unconjugated
	Purification	Antigen Affinity Purified
	Dilution Range	WB 1:1000 ICC 1:50
	Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
	Isotype	IgG
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	3728
	Gene Symbol	JUP
	UniProt ID	PLAK_HUMAN
	Specificity	This antibody was cross-adsorbed to both phospho-tyrosine coupled to agarose and to dephosphorylated Gamma-Catenin (Tyr-550) peptide before affinity purification using phospho-Gamma-Catenin (Tyr-550) peptide (without carrier). The antibody detects an
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.