

## Anti-CTNND1 antibody (central Region) [M354] (STJA0005310)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ICC/IP          |
| Host / Source | Mouse              |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                     |
| Conjugation         | Unconjugated                                                                                   |
| Purification        | Protein A Purified                                                                             |
| Dilution Range      | WB 1:1000<br>ICC 1:100<br>IP 1:50                                                              |
| Formulation         | PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol                                     |
| Isotype             | IgG1                                                                                           |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                  |                                                                                                                                                       |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID          | <a href="#">1500</a>                                                                                                                                  |
| Gene Symbol      | <a href="#">CTNND1</a>                                                                                                                                |
| UniProt ID       | <a href="#">CTND1_HUMAN</a>                                                                                                                           |
| Immunogen Region | central Region                                                                                                                                        |
| Specificity      | The antibody detects a 110 kDa* protein corresponding to the molecular mass of Delta 1-Catenin on SDS-PAGE immunoblots of human A431 and HUVEC cells. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|