

Anti-CD55 antibody (Extracellular region) [M033] (STJA0005324)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:1000 ICC 1:100 IP 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1604
Gene Symbol	CD55
UniProt ID	DAF_HUMAN
Immunogen Region	Extracellular region
Specificity	Clone M033 was purified using Protein G chromatography. The antibody. The antibody also detects a native ^h human recombinant CD55 protein that includes the extracellular region. The antibody can be used for native western blot, immunoprecipitation, pr

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---