

Anti-CD244 antibody (Extracellular region) [M053] (STJA0005334)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IP/ICC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:500 ICC 1:200 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51744
Gene Symbol	CD244
UniProt ID	CD244_HUMAN
Immunogen Region	Extracellular region
Specificity	Clone M053 mouse monoclonal antibody detects a 110 kDa* protein on SDS-PAGE Native* or denatured immunoblots of human THP1 monocytes. The antibody prefers the native CD244 protein, but also detects denatured CD244. The antibody works for western blot

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---