

Anti-CD47 antibody (Extracellular region) [M054] (STJA0005336)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB-Nat/ICC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:500 ICC 1:200 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	961
Gene Symbol	CD47
UniProt ID	CD47_HUMAN
Immunogen Region	Extracellular region
Specificity	Clone M054 mouse monoclonal antibody detects native human CD47, as well as aldehyde-fixed CD47. The antibody works for staining in immunocytochemistry, capture and detection in ELISA, native immunoprecipitation and western blot. The antibody detects

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---