

Anti-CDC42 antibody [M152] (STJA0005338)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:500 ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64465
Gene Symbol	Cdc42
UniProt ID	CDC42_RAT
Specificity	The antibody detects a 21 kDa* protein in human Jurkat cells and mouse brain. It does not recognize a human RhoA GST fusion protein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---