

## Anti-CLASP1 antibody (C-Term) [KT66] (STJA0005343)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rat
Reactivity	Chicken/Human/Mouse/Rat/Xenopus

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:500 ICC 1:100
Formulation	PBS + 0.05% NaN <sub>3</sub>
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">76707</a>
Gene Symbol	<a href="#">Clasp1</a>
UniProt ID	<a href="#">CLAP1_MOUSE</a>
Immunogen Region	C-Term
Specificity	The antibody detects a 161 kDa* protein corresponding to the molecular mass of CLASP1 on SDS-PAGE immunoblots of human HeLa, mouse brain, and rat PC12 cells. The antibody also detects CLASP1 by immunocytochemistry in HeLa cells.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------