

Anti-Phospho-CORO1B-Ser-2 antibody (STJA0005349)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57175
Gene Symbol	CORO1B
UniProt ID	COR1B_HUMAN
Specificity	The antibody detects a 60 kDa* band corresponding to Coronin-1B on SDS-PAGE immunoblots of human A431 or HeLa cells treated with Calyculin A. This band is not observed in control cells and the band observed in Calyculin A-treated cells can be removed

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---