

Anti-Phospho-CREB1-Ser-133 antibody [M415] (STJA0005351)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/ICC
	Host / Source	Mouse
	Reactivity	Chicken/Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Protein A Purified
	Dilution Range	WB 1:1000 ICC 1:100
	Formulation	PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol
	Isotype	IgG1
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	1385
	Gene Symbol	CREB1
	UniProt ID	CREB1_HUMAN
	Specificity	The antibody detects a 43 kDa* protein corresponding to the molecular mass of CREB on SDS-PAGE immunoblots of human A431, Jurkat, and rat PC12 cells treated with Calyculin A, as well as A431 treated with EGF and PC12 treated with NGF or Forskolin. Th
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.