

## Anti-DPYSL2 antibody (C-Term) (STJA0005358)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ICC/IP          |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                     |
| Conjugation         | Unconjugated                                                                                   |
| Purification        | Antigen Affinity Purified                                                                      |
| Dilution Range      | WB 1:1000<br>ICC 1:50<br>IP 1:100                                                              |
| Formulation         | PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol                                     |
| Isotype             | IgG                                                                                            |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                  |                                                                                                                                                                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID          | <a href="#">1808</a>                                                                                                                                                                                                                                      |
| Gene Symbol      | <a href="#">DPYSL2</a>                                                                                                                                                                                                                                    |
| UniProt ID       | <a href="#">DPYL2_HUMAN</a>                                                                                                                                                                                                                               |
| Immunogen Region | C-Term                                                                                                                                                                                                                                                    |
| Specificity      | The antibody detects a 70 kDa* protein corresponding to the molecular mass of CRMP2 on SDS-PAGE immunoblots of human SKN-SH cells, mouse brain, and rat PC12 cells. This band on western blot can be specifically blocked using CRMP2 (C-terminal region) |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|