

## Anti-Phospho-DPYSL2-Thr-555 antibody (STJA0005360)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1808</a>
Gene Symbol	<a href="#">DPYSL2</a>
UniProt ID	<a href="#">DPYL2_HUMAN</a>
Specificity	This antibody was cross-adsorbed to unphosphorylated CRMP2 (Thr-555) before affinity purification using phospho-CRMP2 (Thr-555) peptide (without carrier). The antibody detects a 70 kDa* protein corresponding to CRMP2 on immunoblots of rat PC12 or mou

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------