

Anti-PPIB antibody [M019] (STJA0005365)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:1000 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5479
Gene Symbol	PPIB
UniProt ID	PPIB_HUMAN
Specificity	Clone M019 detects a 20 kDa* band corresponding to the molecular mass of cyclophilin B on SDS-PAGE immunoblots of human LNCaP, MeWo, MCF7, and MDA-MB-231 cell lysates, as well as full length recombinant cyclophilin B, but does not detect recombinant

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---