

Anti-DIAPH2 antibody (C-Term) (STJA0005368)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81624
Gene Symbol	DIAPH3
UniProt ID	DIAP3_HUMAN
Immunogen Region	C-Term
Specificity	The antibody detects a 135 kDa* protein corresponding to the apparent molecular mass of mDia2 on SDS-PAGE immunoblots of various cells and tissues, including human A431 and HeLa cells, mouse brain tissue and C2C12 cells, rat PC12 cells, and rabbit sp

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---