

## Anti-Phospho-EEF2K-Ser-398 antibody (STJA0005392)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                     |
| Conjugation         | Unconjugated                                                                                   |
| Purification        | Antigen Affinity Purified                                                                      |
| Dilution Range      | WB 1:1000                                                                                      |
| Formulation         | PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol                                     |
| Isotype             | IgG                                                                                            |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|             |                                                                                                                                                                                                                                                            |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID     | <a href="#">29904</a>                                                                                                                                                                                                                                      |
| Gene Symbol | <a href="#">EEF2K</a>                                                                                                                                                                                                                                      |
| UniProt ID  | <a href="#">EF2K_HUMAN</a>                                                                                                                                                                                                                                 |
| Specificity | The antibody was cross adsorbed to unphosphorylated eEF2K (Ser-398) peptide before affinity purification using phospho-eEF2K (Ser-398) peptide (without carrier). This antibody detects a 105 kDa* protein corresponding to the molecular mass of eEF2K on |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|