

## Anti-Phospho-EEF2K-Ser-445 antibody (STJA0005393)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Xenopus

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">29904</a>
Gene Symbol	<a href="#">EEF2K</a>
UniProt ID	<a href="#">EF2K_HUMAN</a>
Specificity	The antibody was cross-absorbed to unphosphorylated eEF2K (Ser-445) before affinity purification using phospho-eEF2K (Ser-445) peptide (without carrier). This antibody detects a 105 kDa* protein corresponding to the molecular mass of eEF2K on SDS-PAG

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---