

Anti-Phospho-EGFR-Ser-1142 antibody (STJA0005396)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	1956
Gene Symbol	EGFR
UniProt ID	EGFR_HUMAN
Specificity	This antibody was cross-adsorbed to unphosphorylated EGFR (Ser-1142) coupled to agarose then affinity purified using phospho-EGFR (Ser-1142) peptide (without carrier). The antibody detects a 180 kDa* protein on SDS-PAGE immunoblots of human A431 cell
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.