

## Anti-EPHA2 antibody (Extracellular region) [M060] (STJA0005402)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB-Nat/ICC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1969</a>
Gene Symbol	<a href="#">EPHA2</a>
UniProt ID	<a href="#">EPHA2_HUMAN</a>
Immunogen Region	Extracellular region
Specificity	Clone M060 mouse monoclonal antibody detects a 125kDa* protein on SDS-PAGE Native" immunoblots of human NCI-H2052 and MDA-MB-231 carcinomas, as well as the EphA2 C-terminal 60 kDa fragment in human A431 and A549. The antibody detects native EphA2, a

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---