

Anti-Phospho-EPHA4-Tyr-602 antibody (STJA0005406)

GENERAL INFORMATION

| | |
|---------------|-------------------------|
| Product Type | Primary antibodies |
| Applications | WB/IHC |
| Host / Source | Rabbit |
| Reactivity | Chicken/Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Polyclonal |
| Conjugation | Unconjugated |
| Purification | Antigen Affinity Purified |
| Dilution Range | WB 1:1000 IHC 1:100 |
| Formulation | PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|--|
| Gene ID | 2043 |
| Gene Symbol | EPHA4 |
| UniProt ID | EPHA4_HUMAN |
| Specificity | This antibody was cross-adsorbed to unrelated phosphotyrosine peptides before affinity purification using EphA4 (Tyr-602) peptide. The purified antibody detects a 120 kDa* band corresponding to EphA4 in Western blots of HUVEC cells treated with perva |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|