

Anti-CCDC88A antibody (C-Term) [M012] (STJA0005427)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IHC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:250 IHC 1:50 ICC 1:50 IP 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55704
Gene Symbol	CCDC88A
UniProt ID	GRDN_HUMAN
Immunogen Region	C-Term
Specificity	Clone M012 mouse monoclonal antibody detects a 250 kDa* protein on SDS-PAGE immunoblots of human A431 cells and human brain tissue. The antibody also works for immunoprecipitation, immunohistochemistry, and immunocytochemistry.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---