

Anti-GSK3B antibody (N-Term) [M131] (STJA0005433)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/ICC/IHC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500 IHC 1:250 ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	2932
Gene Symbol	GSK3B
UniProt ID	GSK3B_HUMAN
Immunogen Region	N-Term
Specificity	The antibody detects a 46 kDa* protein corresponding to the apparent molecular mass of GSK-3 Beta on SDS-PAGE immunoblots of rabbit fibroblasts, mouse brain, as well as human SKN-SH and A431 cells.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.