

## Anti-N/A antibody (C-Term) [M050] (STJA0005435)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IHC/IP/FC
Host / Source	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:5000 IHC 1:100 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Immunogen Region	C-Term
Specificity	Clone M050 detects recombinant proteins with C-terminal His-tag epitope. The antibody binds 6 and 10 His-tag epitopes on the C-terminal side of a protein, but does not bind N-terminal His-tag epitopes. The antibody works for western blot, immunopreci

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---