

Anti-HSPD1 antibody (N-Term) [M438] (STJA0005440)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3329
Gene Symbol	HSPD1
UniProt ID	CH60_HUMAN
Immunogen Region	N-Term
Specificity	This antibody detects a 60 kDa* protein corresponding to the apparent molecular mass of Hsp60 on SDS-PAGE immunoblots of human Jurkat and rat A7r5 cells. In immunocytochemistry, anti-Hsp60 specifically stains mitochondria in paraformaldehyde fixed an

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---