

## Anti-ITGB4 antibody (Cytoplasmic) [M126] (STJA0005451)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:250
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3691</a>
Gene Symbol	<a href="#">ITGB4</a>
UniProt ID	<a href="#">ITB4_HUMAN</a>
Immunogen Region	Cytoplasmic
Specificity	This antibody detects a 200kDa* protein corresponding to the molecular mass of Integrin Beta 4 on SDS-PAGE immunoblots of human A431 cells.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---