

Anti-MAPK8 antibody (C-Term) [M267] (STJA0005455)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IHC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 IHC 1:100 ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5599
Gene Symbol	MAPK8
UniProt ID	MK08_HUMAN
Immunogen Region	C-Term
Specificity	This antibody detects a 46 kDa* protein corresponding to the apparent molecular mass of JNK1 on SDS-PAGE immunoblots of human A431 and HeLa, as well as rat PC12 cells.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---