

Anti-IMMT antibody [M027] (STJA0005467)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10989
Gene Symbol	IMMT
UniProt ID	MIC60_HUMAN
Specificity	Clone M027 detects a 90 kDa* band corresponding to the molecular mass of mitofilin on SDS-PAGE immunoblots of native or denatured human PC3, MDA-MB-231, A431, A549, and MeWo cell lysates, as well as human breast, skin, and brain tissues. In addition,

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---