

Anti-NRP1 antibody (a1 CUB Domain) [M534] (STJA0005483)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IHC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:250 IHC 1:100 ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8829
Gene Symbol	NRP1
UniProt ID	NRP1_HUMAN
Immunogen Region	a1 CUB Domain
Specificity	The antibody detects a 130 kDa* protein corresponding to the apparent molecular mass of neuropilin-1 on SDS-PAGE immunoblots of human recombinant neuropilin-1 and in human PC3 cells.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---