

Anti-CD274 antibody (Extracellular region) [M080] (STJA0005524)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/ICC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500 ICC 1:200 IP 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
UniProt ID	PD1L1_HUMAN
Immunogen Region	Extracellular region
Specificity	Clone M080 mouse monoclonal antibody detects a 45-55 kDa* protein on SDS-PAGE Native* or denatured immunoblots of human recombinant PD-L1 extracellular region, and human lung tissue. The antibody works for western blot, immunoprecipitation, immunocytochemistry.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---