

Anti-PRKCA antibody [M499] (STJA0005532)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/ICC
	Host / Source	Mouse
	Reactivity	Chicken/Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Antigen Affinity Purified
	Dilution Range	WB 1:1000 ICC 1:100
	Formulation	PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol
	Isotype	IgG2a
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	5582
	Gene Symbol	PRKCG
	UniProt ID	KPCG_HUMAN
	Specificity	This antibody detects 80-82kDa* proteins corresponding to the molecular mass of PKC Alpha, PKC Beta, and PKC Gamma on SDS-PAGE immunoblots of mouse brain lysates. Similar results were observed in human and mouse lysates. This antibody detects recombi
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.