

Anti-Phospho-PRKCA-S657/Y658 antibody (STJA0005533)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IP
Host / Source	Rabbit
Reactivity	Human/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5578
Gene Symbol	PRKCA
UniProt ID	KPCA_HUMAN
Specificity	This antibody detects an 82kDa* protein corresponding to the molecular mass of phosphorylated PKC Alpha on SDS-PAGE immunoblots of neonatal rat brain lysate. Similar results were observed in human SKN-SH, endothelial, and HeLa cells, as well as rabbi

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---